

SNPsig® EscapePLEX™ SARS-CoV-2

Are you looking for a multi-assay to identify the aggressive variants of SARS-CoV-2 affecting your local communities?

Research
Use Only

Look no further, with EscapePLEX™ the management of SARS-CoV-2 just got **BETTER** and **CLEARER**.

DETECTION PROFILE

- The E484K, K417N, K417T and P681R mutations which are found in Beta, Gamma, Delta and Delta Plus Variants of Concern
- Plus ORF1ab and M gene targets for confirmation of positive SARS-CoV-2

KEY FEATURES

- Use as a reflex test with positive SARS CoV-2 samples.
- Complementary to NGS.
- A faster and more cost-effective alternative to ID E484K, K417N, K417T and P681R mutations.
- Streamline your workflow for easy use with readily prepared Master Mix.
- Open platform instrumentation.
- Shipping at ambient temperature.

ORDERING INFORMATION

PRODUCT CODE: R00152
PACK SIZE: 96 reactions
(Master Mix included)

Primerdesign Ltd.
York house, School Lane, Chandler's Ford S053 4DG
United Kingdom

www.primerdesign.co.uk
www.novacyt.com

Primerdesign is part of the Novacyt Group

Contact us today to place your order:
Telephone: +1 786-503-5180
enquiries@primerdesign.co.uk

